

# Upgrade of HIT's slow extraction system and the first experience in medical operation

*Monday, 12 February 2024 15:20 (20 minutes)*

The Heidelberg Ion-Beam Therapy Centre HIT uses the RF-KO slow extraction method to extract the particles from the synchrotron. To improve the spill quality of the extracted beam a new RF-signal was investigated which increases the R-value from 92.5% to 97.5%. The signal is a multiband RF signal broadened with a random BPSK at 3 frequency bands. The new KO-DDS which generates the signal is in medical use since October 2023 and the first experience are presented in this talk.

**Primary author:** FELDMEIERS, Eike (HIT)

**Co-authors:** PETERS, Andreas (Heidelberger Ionenstrahl-Therapiezentrum (HIT)); SCHÖMERS, Christian (Heidelberg Ion Beam Therapy Centre); CORTÉS GARCÍA, Edgar Cristopher (DESY); HUN, Marcel (Heidelberg Ion-Beam Therapy Center); HABERER, Thomas (Heidelberg Ionbeam Therapy Center)

**Presenter:** FELDMEIERS, Eike (HIT)

**Session Classification:** Spill Ripples & Beam Quality