

Electro-Magnetic Pulse Measurements at the VEGA Laser Facility

M. Ehret, IMPULSE Project at CLPU, Villamayor, Spain


J. I. Apiñaniz, A. Curcio, J. Hernández, J. A. Pérez, P. Puyuelo Valdés and G. Gatti, CLPU, Villamayor, Spain

Goals


- ★ systematic record of the EMP strength inside the interaction chamber and in Target Area (TA) for gas/solid/liquid targets driven by 100ds of TW to PW relativistic laser pulses
- ★ correlation of set-up and EMP spectrum
- ★ user support and failure analysis

Set-up and Methodology

VEGA-2 Experiment Station



VEGA-3 Experiment Station




- ★ E- and B-field antenna signal recorded on GHz oscilloscope (calibrated circuit!)
- ★ analytical calculation of modes in cylindrical and rectangular cavities
- ★ towards high repetition rate online analysis with python software 'pyWA'


Gas Target

- ★ driver laser VEGA-2 at (3–5) J in 30 fs to several 10^{19} W/cm² by F/13
- ★ 8 mm FWHM gas jet target at 10^{18} /cc He


Cylindrical Cavity Modes (B-Field Antenna)




Potential Difference of Ground and Table S



External Spectrum (E-Field Antenna at A)



Typical EMP in TA (E-Field Antenna)





→ 100s of μ T B-field amplitude (≈ 10 stronger than Earth's magnetic field), and 100s of V/m E-field (maximum fields sent by GSM mobile phones)

Experimental Results


Solid Target

- ★ driver laser VEGA-3 at (10–30) J in 30 fs to several 10^{20} W/cm² by F/10
- ★ solid density target of 3 μ m thickness

Cavity Modes of Rectangular VEGA-3 Interaction Chamber (B-Field Antenna at (a,h) = 17.5 cm)



add thick LiF obstacle behind main target



→ spatio-temporal evolution of laser-driven charge distributions has an EMP fingerprint

Perspectives


- ★ careful noise subtraction (mobile phones, WiFi, etc.)
- ★ determination of power in EMP
- ★ identification of applications (outcoupling of modes, etc.)
- ★ measure & log the EMP in TA with an antenna array
- ★ identify 'calm' bands for experiment specific RF measurements with sensitive B-dot and D-dot probes

Toolbox

- ★ calibrated commercial solutions

Magnetic Field

- MDF9400 (Aaronia)
- ★ 9 kHz – 400 MHz



Electric Field

- OmniLOG 30800 (Aaronia)
- ★ 300 MHz – 8 GHz

