

II. Physikalisches
Institut



JUSTUS-LIEBIG-
UNIVERSITÄT
GIESSEN



HELMHOLTZ
GEMEINSCHAFT



FRS Ion Catcher: First Results and Future Perspectives

Wolfgang Plaß

GSI Helmholtzzentrum für Schwerionenforschung, Darmstadt
II. Physikalisches Institut, Justus-Liebig-Universität Gießen, Germany



Overview

- Motivation
- Cryogenic Stopping Cell for the Super-FRS at FAIR
- Multiple-Reflection Time-of-Flight Mass Spectrometer
- On-line Commissioning at the FRS Ion Catcher at GSI
- Conclusions

Motivation: Low Energy Branch of the Super-FRS

LEB: High-precision experiments with in-flight separated exotic nuclei almost at rest, (production by projectile fragmentation / fission)

- universal and fast production
- high selectivity
- cooled exotic nuclei



Eur. Phys. J. Special Topics 183 (2010) 1

MATS (Precision Measurements of very short-lived nuclei using an Advanced Trapping System for highly charged ions)


- High accuracy mass measurements
- In-trap conversion electron and alpha spectroscopy
- Trap assisted spectroscopy

LaSpec (Laser Spectroscopy)

- Collinear laser spectroscopy of ions and atoms
- β -NMR
- Resonance ionization spectroscopy


Stopping Cell Principles

Range Bunching



H. Weick et al., NIM B 164 (2000) 168

Stopping Cell



- High-energy ions stopped in noble gas
- Stopped ions transported using DC and RF fields to exit-hole
- Extraction by gas flow


Stopping Cell for Relativistic Exotic Nuclei

First Generation Stopping Cell (S258)

Successful proof-of-principle

Suffered from:


1. Poor stopping efficiency
2. Presence of impurities



Second Generation Stopping Cell (S411)

Solutions now implemented:

1. High-density operation
2. Cryogenic operation




Stopping Cell: Cryogenic Operation

New concept: Operate He-filled stopping cell at cryogenic temperature (~ 70 K)

Advantages

- Ultra-pure helium (freezing-out of contaminants)
 - Ideal for ion survival
 - No formation of molecules/adducts
- Reduced radial ion diffusion
- 2+ charge state (?) \rightarrow shorter extraction times
- Reduced requirements for cleanliness \rightarrow easier, more flexible construction
- Operational reliability




Transport efficiency
of α -decay recoil ions
in a closed gas cell

P. Dendooven et al., NIM A 558 (2006) 580

S. Purushothaman et al., NIM B 266 (2008) 4488

Cryogenic Stopping Cell Design




M. Ranjan et al., Europhys. Lett. 96 (2011) 52001
 M. P. Reiter, Master Thesis, Justus-Liebig-Universität Gießen (2011)

Developed in Collaboration

 Kernfysisch Versneller Instituut
 JUSTUS-LIEBIG-UNIVERSITÄT GIESSEN
 university of groningen
 GSI



Cryogenic Stopping Cell: RF Carpet



Off-line Performance Study: Efficiency

S. Purushothaman et al.

Measurement with ^{223}Ra source



Maximum survival and extraction efficiency for ^{219}Rn : ~ 30%

Value is close to the assumed efficiency value for stopping the recoils as ions

(cf. P. Dendooven et al., NIM A 558 (2006) 580; S. Purushothaman et al, NIM B 266 (2008) 4488;

S. Eliseev et al., NIM B 258 (2007) 479.)



Maximum pressure reached: 100 mbar, 90 K

→ 330 mbar room-temperature equivalent; 5 mg/cm² (He)

In the future: RF = 150 V → 1 bar room-temperature equivalent achievable

Off-line Performance Study: Extraction Time

Extraction time measurement with pulsed ^{223}Ra source



Pressure = 50 mbar, temperature = 100 K
DC field = 21 V/cm, distance = 1 m


→ Experiment $T_{\text{extr}} \sim 45$ ms
Theory ($K_0=15$ cm²/Vs) $T_{\text{extr}} \sim 47$ ms

Scaled mean extraction time
(extraction from middle of cell)


$$T_{\text{extr}} \sim 23 \text{ ms}$$

T. Dickel, S. Purushothaman, et al.

Multiple-Reflection Time-of-Flight Mass Spectrometer




Isobar Separation Mode




$m/\Delta m > 10^5$

High Resolution Mode



$m/\Delta m > 10^5$,
 Mass Accuracy 10^{-6} - 10^{-7}

Broadband Mode



Full Mass Range,
 $m/\Delta m \sim 2000$


W.R. Plaß et al., NIM B 266 (2008) 4560

MR-TOF-MS: Mass Resolving Power

Isobar doublet at mass 78 u

- pyridine fragment (C_5H_5N)
- benzene (C_6H_6)

256 Turns, $TOF_{C_5H_4N} = 9.2$ ms, $m/\Delta m = 160,000$



No re-tuning, only delays changed
Data taken within 6 minutes

T. Dickel


MR-TOF-MS: Isobar Separator

Combination of MR-TOF-MS and Bradbury-Nielsen Gate yields powerful isobar separator

- Proposed in 2005 (W.R. Plaß et al.)
- First demonstrated in 2006 W.R. Plaß et al., NIM B 266 (2008) 4560

Isobar separation demonstrated for C_6H_6 and $^{13}C^{12}C_5H_5$

- Intensity ratio 170:1
- $\Delta m = 4$ mu
- up to 10^4 ions/cycle
→ 10^6 ions/s



T. Dickel, PhD Thesis Justus-Liebig-Universität Giessen, 2010

Performance Characteristics of the MR-TOF-MS

Universal mass spectrometer and mass separator
(works for all elements, stable and unstable ions)

Mass Resolving Power

600,000

Mass Range

Broadband: ~% - full

Mass Measurement Accuracy

$\sim 10^{-7}$

Transmission efficiency

up to 70%

Measurement Duration

~ 10 ms

Ion Capacity

$> 10^6$ ions / s

Sensitivity

~ 10 ions

Dynamic Range

$> 10^4$

**World-wide unique combination of performance characteristics,
ideally suited for the LEB of the Super-FRS**

Further performance improvements are underway..


T. Dickel

FRS Ion Catcher

Fragment separator FRS at GSI is the ideal instrument for testing Super-FRS developments

FRS Ion Catcher


- Test bench for the cryogenic stopping cell and the MR-TOF-MS
- Potential for high-precision experiments with stopped projectile and fission fragments
 - Direct mass measurements
 - Mass-selected decay spectroscopy
 -



W.R. Plaß et al., GSI Scientific Report 2010 (2011) p. 137


Challenges

- Limited space at the final focal plane
- Setup cannot be installed permanently
- Short setup times
- Limited amount of beam time



Projectile Fragment Separator FRS

FRS Ion Catcher Experiment in October 2011



Installation at the FRS

Transport of stopping cell and MR-TOF-MS from off-line location to FRS




The components just fit into the available space at the final focal plane of the FRS



Time between MR-TOF-MS spectra at the offline location and the FRS: 4 hours!


First mobile high-performance MS

FRS Ion Catcher Setup



Range Distribution and Stopping Efficiency

Range distribution of ^{223}Th




- Range distribution: $\sigma = 14 \text{ mg/cm}^2$ (Al) $\rightarrow \sigma = 11 \text{ mg/cm}^2$ (He)
- Stopping cell gas density: 100 mbar, 100 K $\rightarrow 0.05 \text{ mg/cm}^3$
Areal weight of stopping cell: 5 mg/cm^2 (He)

Stopping efficiency = $\sim 16\%$ (Preliminary)

Stopping Cell Extraction Efficiency

Si detector spectrum of extracted projectile fragments




First on-line operation of a cryogenic stopping cell

Efficiencies for ^{223}Th (preliminary)


- Overall efficiency ~ 8%
- Survival and extraction efficiency ~ 50%

Measurement with the MR-TOF-MS

Measurement with ^{223}Ra source




Mass Spectrum



α Spectrum

Broadband mass spectrometry is a necessity for quick and reliable operation of a stopping cell




J. Ebert, Master Thesis, Justus-Liebig-Universität Gießen (2011)

Conclusions


Stopping cell for the Super-FRS and the FRS Ion Catcher

- Cryogenic, high density operation
- Commissioned off-line and on-line
- Preliminary performance values:
Stopping efficiency ~ 16%, survival and extraction efficiency ~ 50%
Extraction time ~ 20 ms



High-performance multiple-reflection time-of-flight mass spectrometer

- Direct mass measurements
- High-resolution mass separator
- Diagnostics tool: identification and quantification
- World-wide unique combination of performance characteristics, ideally suited for current and next generation facilities



Outlook

Next steps

- Commission cryo-cooler based cooling system for the stopping cell
- Systematic study of the cryogenic stopping cell (intensity limitations and temperature / cleanliness)
- Mass measurements with the MR-TOF-MS
- Increase stopping efficiency even further (higher densities)



Acknowledgements

FRS Ion Catcher / S411 Collaboration

P. Dendooven, T. Dickel, J. Ebert, A. Estrade, F. Farinon, H. Geissel, E. Haettner, C. Jesch, N. Kalantar-Nayestanaki, R. Knöbel, J. Kurcewicz, J. Lang, I.D. Moore, C. Nociforo, H. Penttilä, S. Pietri, W.R. Plaß, A. Prochazka, S. Purushothaman, M. Ranjan, M.P. Reiter, S. Rinta-Antila, C. Scheidenberger, M. Takechi, H. Weick, J.S. Winfield



C. Jesch

IONAS Group at JLU Gießen

S. Ayet, U. Czok, T. Dickel, M. Diwisch, J. Ebert, H. Geissel, F. Greiner, E. Haettner, C. Horbach, C. Jesch, N. Kuzminchuk, R. Knöbel, J. Lang, M. Petrick, W.R. Plaß, M.P. Reiter, A.-K. Rink, C. Scheidenberger, B. Sun, W. Lippert, B. Wohlfahrt, M. Yavor

IONAS Alumni

A. Becker, L. Chen, Z. Di, S. Eliseev, B. Fabian, T. Fleckenstein, W. Kinsel, F. Lautenschläger, K. Reinheimer, D. Schäfer, T. Schäfer, A. Simon, Z. Wang, J. Werner



**Funding: Univ. Groningen and GSI,
HGF and GSI (VH-NG 33), GSI F&E (GIMET2)
BMBF (06GI185I, 06GI9114I)**

