

RIKEN GARIS for Superheavy Element Chemistry



Nishina Center, RIKEN
Hiromitsu Haba




CONTENTS

1. Introduction
2. Present status of GARIS
 - 2.1. Developments in 2008
 - Rotating ^{248}Cm target and new gas-jet chamber
 - 2.2. Production of ^{261}Rf and ^{265}Sg
3. Chemistry programs
4. Summary and perspectives





1. Introduction

RIKEN GARIS as a pre-separator for SHE chemistry



Startup of the SHE chemistry in RIKEN

Model experiments with a prototype gas-jet transport system



Powerful tool for SHE chemistry

Extremely low background condition

Beam-independent high gas-jet efficiency

For the future SHE chemistry

→ ^{248}Cm -based hot fusion reactions

Z	Reaction	σ (pb)
104	$^{248}\text{Cm}(^{18}\text{O},5n)^{261}\text{Rf}$	13000
105	$^{248}\text{Cm}(^{19}\text{F},5n)^{262}\text{Db}$	1500
106	$^{248}\text{Cm}(^{22}\text{Ne},5n)^{265}\text{Sg}$	240
107	$^{248}\text{Cm}(^{23}\text{Na},5n)^{266}\text{Bh}$	40*
108	$^{248}\text{Cm}(^{26}\text{Mg},5n)^{269}\text{Hs}$	7

* from the σ(5n) vs. Z systematics

In this presentation

1. Present status of GARIS

Developments in 2008

Rotating ^{248}Cm target

New gas-jet chamber

Chemistry laboratory

→ **Production of ^{261}Rf and ^{265}Sg**

2. Chemistry programs at RIKEN


3. Future plans

2. Present status of GARIS



2.1. Developments in 2008

Rotating ^{248}Cm target

- Purification with ion exchange
- Electrodeposition
0.54 mg of ^{248}Cm in 20 μL of 0.2 M HNO_3 + 5.5 mL 2-propanol
1000 V x 11 mA/cm² for 10 min
→ **280 $\mu\text{g}/\text{cm}^2$ $^{248}\text{Cm}_2\text{O}_3$ on 2.0 μm Ti backing foil**



Water-cooled cell for electrodeposition




New gas-jet chamber

Size of the focal plane: $\Phi 60\text{ mm} \rightarrow \Phi 100\text{ mm}$

Support mesh: $\Phi 4\text{-mm circular holes (72\%)} \rightarrow \Phi 2\text{-mm circular holes (78\%)}$

Mylar window: $1.1\text{ }\mu\text{m} \rightarrow 0.5\text{ }\mu\text{m thickness}$




New gas-jet chamber@GARIS

Chemistry laboratory



Chemistry laboratory
Focal plane of GARIS



2.2. Production of ^{261}Rf and ^{265}Sg

Experimental conditions for ^{261}Rf





	$^{248}\text{Cm}(^{18}\text{O},5n)^{261}\text{Rf}$	$^{\text{nat}}\text{Gd}(^{18}\text{O},xn)^{169}\text{Hf}$
Cross section	13 nb ¹⁾	170 mb ²⁾
Beam energy (MeV)	95.4	←
Beam intensity (pμA)	6	←
Target on 2 μm Ti (μg/cm ²)	250	300
ER recoil energy (MeV)	6.3	9.6
Magnetic rigidity (Tm)	1.58-2.16	1.48-1.63
GARIS He (Pa)	32	←
Mylar window (μm)	0.54	←
Support mesh (%)	78	←
Gas-jet He (kPa)	49	←
He flow rate (L/min)	2	←
KCl generator (°C)	620	←

1) Nagame et al.: J. Nucl. Radiochem. Sci. **3**, 85 (2002).



2) Calculated with PACE4.

Experimental setup


GARIS





$^{248}\text{Cm}(^{18}\text{O},5n)^{261}\text{Rf}$



No background peaks at > 8 MeV
161 α (58 α - α) / 8.2-h meas.
 $B_p = 1.75 \pm 0.02$ Tm, $\Delta B_p/B_p = 12.7 \pm 1.9\%$
GARIS eff.: $7.8 \pm 1.7\%$ for $\Phi 100$ mm ($\sigma = 13$ nb)
Gas-jet eff.: $52 \pm 12\%$



$^{nat}\text{Gd}(^{18}\text{O},xn)^{169}\text{Hf}$



Mixed $^{248}\text{Cm}/^{nat}\text{Gd}$ target

→ Simultaneous chemical experiments with ^{261}Rf and ^{169}Hf

$B\rho = 1.57 \pm 0.01$ Tm for ^{169}Hf

Change of the magnet settings of GARIS: ~1 min for $^{261}\text{Rf} \leftrightarrow ^{169}\text{Hf}$

Production of ^{265}Sg with the GARIS/gas-jet system


Oct. 1–6, 2008

Experimental conditions for ^{265}Sg

Reaction	$^{248}\text{Cm}(^{22}\text{Ne},5n)^{265}\text{Sg}$
Cross section	200-300 pb [*]
Beam energy (MeV)	117.9
Beam intensity (μA)	3
Target on 2 μm Ti ($\mu\text{g}/\text{cm}^2$)	250
Recoil energy (MeV)	9.4
Magnetic rigidity (Tm)	1.73, 1.94, 2.05, 2.16
GARIS He (Pa)	32
Mylar window (μm)	0.65
Support grid (%)	78
Gas-jet He (kPa)	48
He flow rate (L/min)	2
KCl generator ($^\circ\text{C}$)	600

*Düllmann and Tüler: Phys. Rev. C 77, 064320 (2008).

$^{248}\text{Cm}(^{22}\text{Ne},5n)^{265}\text{Sg}$



14 correlations (35 α /fission events) on ^{265}Sg , ^{261}Rf , and ^{257}No
 $B\rho = 2.07 \pm 0.01 \text{ Tm}$, $\Delta B\rho/B\rho = 8.4 \pm 1.1\%$

3. Chemistry programs

Workshop on SHE Chemistry at RIKEN (Dec. 11–12, 2007, KUR)

What chemistries should be studied at RIKEN?

GARIS@RILAC

1. Decay studies of Rf and Sg isotopes for chemical investigations
(H. Kudo of Niigata Univ.)
2. Solvent extraction of Sg with micro-chemical chip and LS
(A. Shinohara of Osaka Univ.)
3. Electrochemistry of Sg and Bh with flow electrolytic cell
(A. Toyoshima of JAEA)

@AVF Cyclotron

1. Electrochemistry of Md with flow electrolytic cell (A. Toyoshima of JAEA)
2. Decay studies of Sg isotopes using gas-phase chemistry
(T. Sato of JAEA).
3. Reversed-phase extraction chromatography of Rf with AIDA
(A. Yokoyama of Kanazawa Univ.)

Workshop on SHE Chemistry@RIKEN in 2008 (Nov. 12, 2008)

4. Summary and perspectives

Present status of RIKEN GARIS

- Developments in 2008
 - Rotating ^{248}Cm target
 - New gas-jet chamber
 - Chemistry laboratory
- Production of SHEs for chemistry experiments
 - $^{248}\text{Cm}(^{18}\text{O},5n)^{261}\text{Rf}$ and $^{248}\text{Cm}(^{22}\text{Ne},5n)^{265}\text{Sg}$

Future plans

- New separator GARIS II for hot fusion reactions
 - $Q_vD(30^\circ)Q_hQ_vD(7^\circ)$, 20-msr solid angle, $B\rho_{\max} = 2.3 \text{ Tm}$, 5.1-m length
 - Installation in 2009 and commissioning from 2010
- Chemistry experiments
 - Workshop on SHE Chemistry@RIKEN in 2008 (Nov. 12, 2008)
 - The 4th PAC Meeting at RIBF (Jan., 2009)

Collaborators

RIKEN

T. Ichikawa, D. Kaji, Y. Kudou, K. Morimoto, K. Morita, K. Ozeki, N. Sato, A. Yoneda, and A. Yoshida

JAEA

M. Asai, Y. Kasamatsu, Y. Nagame, T. Sato, A. Toyoshima, and K. Tsukada

Osaka Univ.

H. Fujisawa, H. Kikunaga, Y. Komori, K. Ooe, A. Shinohara, W. Yahagi, and T. Yoshimura

Niigata Univ.

S. Goto, T. Kawasaki, and H. Kudo

Kanazawa Univ.

M. Araki, T. Nanri, and A. Yokoyama

Thank you!