

Delineation and Target concepts


Waisse WAISSI MD. PhD.

Département de Radiothérapie

Centre Léon Bérard

Génétique Epigénétique et Biologie des sarcomes

Centre de Recherche en Cancérologie de Lyon



CRCL CENTRE DE
RECHERCHE EN
CANCÉROLOGIE
DE LYON



Disclosures

- None

Objectives

- Understand the concept of target volumes in radiation oncology.
- Learn the definitions: **GTV, CTV, PTV, and OAR.**
- Explore the process of contouring and its clinical importance.


Introduction-Why accurate delineation matters ?

- **Precision in Radiotherapy** : Accurate delineation ensures the tumor receives the prescribed dose while sparing healthy tissues to minimize side effects.
- **Consequences of Errors** : Errors in delineation can cause underdosing of tumors reducing efficacy or overdosing normal tissues increasing toxicity.
- **Standards and Protocols** : ICRU guidelines and advanced imaging technologies standardize delineation to ensure treatment consistency and safety.
- **Hierarchical Volume Definition** : Delineation involves defining volumes accounting for visible disease, microscopic spread, motion, and setup uncertainties.

Historical definitions



Target Volume_ Overview



Nicholas G Zaorsky, MD

GTV: gross tumor volume

- All tumor that can be seen on a scan
- For reference, a 1 cm tumor has about 10^9 billion cells

CTV: clinical target volume

- The volume that has to be treated to achieve cure or palliation
- Should be defined before choosing a treatment modality
- Includes "crablike" extensions of the cancer, which contain:
 - microscopic cancer cells that cannot be seen on a scan
 - "Elective" lymph nodes that we think the cancer spread to
- $CTV = GTV + \text{a margin of where we think microscopic disease is located}$

IM: internal margin

- Variations in size, shape, and position of CTV relative to anatomic reference points (e.g. from breathing)

ITV: internal target volume

- Volume that encompasses movement of the CTV and the IM
- $ITV = CTV + IM$

SM: setup margin

- Uncertainties in patient positioning and alignment of therapeutic beams during treatment planning, and all treatment sessions


PTV: planning target volume

- Volume that encompasses movement of the CTV, the IM, and errors in setup / position of the patient and penumbra
- $PTV = CTV + IM + SM$

Historical Evolution of volume definitions

- **Initial Volume Standardization** : The 1978 ICRU report first distinguished between target and treated volumes, initiating volume standardization in radiotherapy.
- **Introduction of GTV, CTV, PTV** : ICRU Report 50 (1993) introduced Gross Tumor Volume, Clinical Target Volume, and Planning Target Volume for clearer dose specifications.
- **Refinement with ITV and PRV** : ICRU Report 62 (1999) added Internal Target Volume and Planning Organ at Risk Volume to address motion and conformal techniques.
- **Adaptation for IMRT** : ICRU Report 83 (2010) adapted volume definitions for intensity-modulated radiation therapy, recommending volume-based dose prescriptions

Example



GTV – Gross Tumor Volume

- **Definition of GTV** : GTV includes all visible tumor tissue and involved lymph : nodes identified via clinical exams and imaging.
- **Importance of Accurate Delineation** : Precise GTV delineation impacts the size and shape of Clinical and Planning Target Volumes, essential for effective treatment.
- **Challenges in GTV Definition** : Despite imaging advances, inter-observer variability poses challenges in consistently defining GTV boundaries.
- **Planning Without Visible Tumor** : When no tumor is visible post-surgery, treatment planning starts with **Clinical Target Volume** instead of GTV.


CTV – Clinical Target Volume

- **Definition of CTV** : CTV includes the gross tumor volume and surrounding areas at risk for microscopic disease spread.
- **CTV Margin Determination** : CTV margins are determined using clinical judgment based on pathological data and tumor failure patterns.
- **Margin Size Variation** : CTV margins vary by tumor type and location, generally from 5 mm to 20 mm, depending on cancer type.
- **Importance of Accurate Delineation** Accurate CTV delineation balances tumor control with minimizing damage to healthy tissue


Volume definition

CTV T = CTV Chest wall R

CTV N = nL4+nL2+IP+IMC



Volume definition CTV N




Gregoire et al. 2013

ITV-Internal Target Volume

- **Purpose of ITV** : ITV accounts for physiological changes in tumor shape and position during treatment, ensuring consistent coverage.
- **Physiological Variations** : Respiratory motion, cardiac activity, and organ filling cause tumor displacement up to 15 mm, influencing ITV design.
- **Imaging Techniques** : Four-dimensional CT (4DCT) captures tumor motion over time to accurately define the ITV for treatment planning.

ITV




BTV-Biological target Volume

- **Concept of BTV** : BTV integrates functional and metabolic imaging like PET to delineate biologically active tumor regions for targeted therapy.
- **BTV Definition Methods** : Methods include fixed SUV thresholds, percentage of SUVmax, and advanced algorithms such as PERCIST for defining BTV.
- **Challenges in BTV** : Challenges include variability in techniques and spatial mismatches with conventional tumor volumes affecting accuracy.
- **Benefits of PET Integration** : Integrating PET improves delineation accuracy and reduces inter-observer variability, supporting personalized treatment.

BTV

BIOLOGICAL TARGET VOLUME?



- PET
- F-miso
- Hypoxia



- MRI/MRS
- Choline/citrate
- Tumor burden



- PET
- IUDR
- Tumor growth



Biological
eye view



PTV – Planning Target Volume

- Definitions and purpose :
 - PTV ensures prescribed dose coverage of the CTV despite setup and treatment uncertainties.
- Margin Calculation Formula
 - The van Herk formula calculates margins using systematic and random error components for dose accuracy.
- Margin Variability factors
 - PTV margins depend on anatomical sites and advanced imaging for smaller, precise treatment margins.

Conclusions : references tables

VOLUME TYPE	DEFINITION	KEY CONSIDERATIONS
GTV	Visible tumor on imaging	No margin; direct imaging
CTV	GTV + microscopic disease margin	5–20 mm based on tumor type
ITV	CTV + internal margin for motion	Accounts for respiration and organ filling
PTV	CTV/ITV + setup and uncertainty margins	Van Herk formula: $M = 2.5\Sigma + 0.7\sigma$
BTV	Functional/metabolic imaging-based volume	PET SUV thresholds or advanced algorithms
TV	Volume receiving prescribed dose	Defined by 95% isodose line

Target Volume Concepts – Overview

