


































Particle therapy machines

Christian Graeff c.graeff@gsi.de

Project funded by



Federal Department of Economic Affairs, Education and Research EAER State Secretariat for Education, Research and Innovation SERI





Agenda




- Fundamentals of particle beam production
- Accelerators: Cyclotrons and Synchrotrons
- Beam delivery: active scanning
- Gantries
- Facility examples

Governing forces in ring accelerators



- The Lorentz force: $\overrightarrow{F} = q(\overrightarrow{E} + \overrightarrow{v} \times \overrightarrow{B})$
- Electrical component $\overrightarrow{F} = q \ \overrightarrow{E}$
- Force acts on object with a charge q
 - in the direction of the electric field E
 - acts on objects at rest (v = 0)
 - accelerates the object
 - independent of the object's mass m
- Objects traversing an electric field will gain kinetic Energy $E_{kin}=qU$


force \vec{F} charge q electric field \vec{E} voltage $\vec{\pmb{U}}$ velocity $\vec{\pmb{v}}$ magnetic field \vec{B}



Governing forces in ring accelerators



- The Lorentz force: $\overrightarrow{F} = q(\overrightarrow{E} + \overrightarrow{v} \times \overrightarrow{B})$
- Magnetic component $\vec{F} = q(\vec{v} \times \vec{B})$:
- ullet acts on a moving object (v>0) with a charge ${f q}$
 - perpendicular to both \overrightarrow{v} and \overrightarrow{B}
 - does <u>not</u> change the magnitude of \vec{v} , but its direction
 - independent of the object's mass m
 - counteracted by centripetal force $\vec{F}=m\cdot \frac{v^2}{R}$, object turns on a circular track



Governing forces in ring accelerators



 Equality of magnetic Lorentz and centripetal force if B and v are perpendicular:

$$q \cdot v \cdot B = m \cdot \frac{v^2}{R}$$


allows to derive <u>magnetic rigidity Brho</u>


$$BR = B\rho = \frac{m}{q}v$$
 [magnetic rigidity in Tm]

accelerator properties:

field strength, bending radius particle beam properties:

mass, charge, speed





Application of Lorentz forces




- Electrical component <u>acceleration</u>
 - can change magnitude of velocity
 - does not scale with velocity
- Magnetic component <u>steering and focusing</u>
 - changes vector, but not magnitude of velocity
 - scales with velocity much more effective at high speeds


Dipoles and Quadrupoles



• Basic forms of magnets in acceleration

Forces on a positive charged ion traveling out of the screen






Focusing

Particles need bigger accelerators...



Protons and ions need much higher energies for sufficient penetration




Magnetic rigidity of this proton beam: 2.2 Tm, of this carbon ion beam: 6.7 Tm

Ring accelerator principle



- To achieve high energies, re-use accelerating cavity
- Requires
 - synchronization of particle rotation and cavity HF
 - some form of magnetic fields to return particles to cavity again and again
- Realizations for particle therapy: cyclotrons and synchrotrons





Nozzle - everybody scans...

Project funded by



Schweizerische Eidgenossenschaft Confédération suisse Confederazione Svizzera Confederaziun svizra



Education and Research EAER

Research and Innovation SERI






Scanned particle beams


- protons and ions have a charge, so can be easily steered with magnets
- Targeting is actually much easier than with photons

- For precise delivery, tumor is divided into slices of equal beam energy
- a thin beam is scanned over the slice to fill it with dose






GSI's Cave M: the first fully active scanning device



IBA scanning nozzle

- Mounted at gantry exit
- short scan distance for protons
- optimized for PBS
 - quadrupoles, vacuum
- strip IC for intensity and position control



Scanning strategy characteristics




Spot scanning



- 0 0 0
- 00000
 - \circ

Raster scanning

- planned spot positions
- delivered spots



- Beam turned on/off for each spot
 - requires fast, frequent beam off technically challenging for synchrotron
 - scanning magnets need to reach position and ,settle' before beam on – fixed minimum spot duration (some ms)
 - synchro-cyclotron: 3+ pulses per spot

- Beam on for transition between spots
 - dose is deposited between spots
 - magnet scan speed and particle rate become highly relevant
 - minimum spot weight to deliver ,most' of the dose on planned position
 - higher min weight in a plan allows for higher particle rate / faster delivery

4.1 beam formation 14



Cyclotrons

Project funded by



Schweizerische Eidgenossenschaft Confédération suisse Confederazione Svizzera Confederaziun svizra









Cyclotron principle







Cyclotron resonance condition



Particles bending radius adheres to Lorentz and centripetal force:

$$q \cdot v \cdot B = m \cdot \frac{v^2}{r} \rightarrow q \cdot B = m \cdot \frac{v}{r} = m \cdot \omega$$

$$\omega_c = \frac{q \cdot B}{m} = 2\pi \cdot f_c$$

- For resonance to the accelerating HF with $f_{HF}=f_c$, flight-time
 - has to be constant and energy-independent
 - equal to half the HF period on each half-turn
- For constant HF, B, this condition can be met for constant m/q
- Not possible for relativistic particles, as m will change with increasing velocity!

RF: radio-frequency; HF: high-frequency


Cyclotron for relativistic particles




- Two options:
- 1. Synchronize HF to particle velocity: synchro-cyclotron
 - B = const; $f_{HF} = f_{particle} \downarrow$
- 2. Increase B field with increasing radius: isochronous cyclotron
 - B \uparrow ; $f_{HF} = f_{particle} = const$
- In both cases, the factor to modify f or B is γ , the relativistic mass increase
- Both are in use for proton therapy
- synchro-cyclotrons can be more compact, as constant high B-field allows smaller radius
- isochronous cyclotrons can produce a quasi-continuous beam: higher duty cycle


IBA synchro-cyclotron cycle






Isochronous cyclotrons














Isochronous (IBA)

C230 isochronous cyclotron


- □ Diameter: 4.34m
- Weight: 210 tons
- Conventional magnets
- ☐ *Bavg*: 1.74T to 2.2T
- ☐ Dee voltage: 55kV to 150kV peak
- ☐ Rf frequency: 106MHz
- Quasi continuous beam
- Average beam current: 300nA
- Energy: 230MeV (0.6MeV @2σ)



super-conducting synchro-cyclotron (IBA

S2C2 synchrocyclotron

- ☐ Diameter: 2.30m
- Weight: <50 tons</p>
- Superconducting magnets
- □ Bavg: 5.64T to 5.24T
- ☐ Dee voltage: 14kV peak
- ☐ Rf frequency: 90MHz to 60MHz 33%
- Pulsed beam at 1kHz rep rate1 ms
- Average beam current: 150nA
- Energy: 230MeV (2.5MeV @2σ)






super-conducting isochronous (Sumitomo)

SC230 isochronous cyclotron


- ☐ Diameter: 2.8m
- □ Weight: 63t
- Superconducting magnets
- □ *Bavg:* 3.2T to 4T
- ☐ Dee voltage: 50 93 kV
- ☐ Rf frequency: 95.2 MHz
- Quasi continuous beam
- ☐ Average beam current: 1000nA
- Energy: 230MeV (<0.5MeV @2σ)</p>



super-conducting synchro-cyclotron (IBA


S2C2 synchrocyclotron

- ☐ Diameter: 2.30m
- Weight: <50 tons</p>
- Superconducting magnets
- ☐ *Bavg*: 5.64T to 5.24T
- ☐ Dee voltage: 14kV peak
- □ Rf frequency: 90MHz to 60MHz 33%
- Pulsed beam at 1kHz rep rate1 ms
- Average beam current: 150nA
- Energy: 230MeV (2.5MeV @2σ)






Cyclotron facility example: OncoRay Dresden




OncoRay (20)
National Center for Radiation Research in Oncology

IBA cyclotron



Energy selection system



OncoRay (2)X

Rational Center for Radiation Research is Oncology

Problem 3.3.: Man berechne die maximale Dicke des Graphitkeiles.

Energy selection: Varian / PSI





Cyclotron energy is fixed – passive variation is unavoidable

'Energy-selection system'


loss of >95% of particles for lowest energy

 'hot' region – both emissions & long-term activation







Beam transport to the patient



Gantry for patient treatment



Treatment room




One-room facility example: Mevion

Mevion Monarch 250

- Clinical since 2013
 - active scanning since 2018
- Synchro-cylcotron @ 9T
 - $E_{max} = 250 \text{ MeV}, r_{extr} = 0.31 \text{ m}$
 - Mass only 20t
 - mounted on isocentric gantry
- Practically no beam transport
- Beam energy ,selection' directly in front of patient
 - secondary neutrons?


≤3.5mSv/Gy (Chen et al. Med Phys 2013)



C400 Carbon ion cyclotron: Caen, France




- Iso-chronous cyclotron for H₂, He, C
 - all with A/Q=2, but still different frequencies f=75 MHz (12C), 75.6 MHz (H₂)
- 7m super-conducting magnet > 700t
- 3 sources mounted on top: fully stripped He, C
- extreme magnet design with very deep valleys, many tuning options
- assembly on-site, tuning & commissioning



C400 Carbon ion cyclotron: Proton extraction



- Proton molecule accelerated:
 H₂+ with A/Q=2
- Stripping foil inserted to extract at 260 MeV
 - mounted on robotic arm, has to be removed for other ions
- dissociates molecule: A/Q=1, changes trajectories to green, smaller circles
- proton-specific beam line



Summary cyclotrons



- Compact, superconducting machines for protons
 - first carbon cyclotrons treating patients soon!
- Provides a constant energy, passive degradation necessary
- Provides a continuously available beam
- Energy changes requires (short) time

- Majority of particle therapy centers
- Many different facility layouts



Synchrotrons

Project funded by



Schweizerische Eidgenossenschaft Confédération suisse Confederazione Svizzera Confederaziun svizra







Synchrotrons




- Idea: Bending magnets keep particles on a single, periodic track of a fixed length L containing one or more HF cavities
 - Consists of several (many!) magnets
- Bending radius r in each magnet is constant, but momentum increases after cavity

$$r = \frac{\gamma m_0 v}{qB} = const$$

=> B has to increase synchronously to $\gamma m_0 v$

• Also cavity HF has to increase with increasing velocity / rotation frequency $f = \frac{v}{L} = f_{HF}/h$ h: harmonic mode


Largest example: LHC @ CERN with $E_{max} = 7 \text{ TeV}$



Synchrotrons need an injector




- In contrast to cyclotron, ions are not directly injected into the cavity
- Minimum velocity is needed to guide the beato the cavity & keep it on the circular track
- Injection chain consists of
 - an ion source
 - a low-energy accelerator, typically a Linac
- Trade-off between injector & synchrotron
 - high injection energy reduces demands on synchrotron, but injector becomes larger
 - Also accelerator chains of multiple synchrotrons possible




Layout Marburg ion beam therapy center





Synchrotron

■ 430 MeV/u ¹²C




Ion source and injector

- 2 ECR ion sources for p/C
- Linac RFQ + IH
- injection at 7 MeV/u








radiofrequency quadrupole





Linac, IH tank

Example: HIT synchrotron


- Injection at 7 MeV/u
- 6 elements
 - dipole bending magnet (BM)
 - 2 quadrupoles for strong focusing (QF, QD)
- Injection & Extraction septum
 - fast electrostatic kickers
- Septum magnet for extraction
- RF cavity












Energy selection synchrotron




- Direct extraction of any desired energy within specs
- Traditionally energy change only possible from spill-to-spill (2-5s)
- NIRS: Energy change on flattop reduced to ~100 ms
 - Limiting factor: adjustment of large dipoles in high-energy beam transport (as in fast cyclotron systems)
 - Regular treatments can be completed with a single spill of max 2 10¹⁰ carbon ions



First carbon gantry in the world




- mass 670t
- upstream scanning before last bending magnet
- 360° rotation
- ~ 1 mm precision in isocenter



Marburg treatment room






- Fixed beam (45°/90°)
- robotic patient couch
- robotic radiography

One-room synchrotron facility example: Protom

- Proton synchrotron
 - Emax= 330 MeV
 - capable of proton radiography
 - d = 6 m
- Was retrofitted into 2 existing Linac bunkers
- Advantages over cyclotrons:
 - more modular, installation through elevator shafts
 - 'cleaner' beam (no degrader)



Summary synchrotrons



- Larger than cyclotrons, but individual magnets much smaller
 - ramping magnets makes superconductivity challenging!
- Desired energy can be directly extracted
- Distinct rhythm of acceleration and extraction
 - multi-energy extraction for much higher duty-cycle

- All currently treating carbon ion facilities
- Some proton facilities

Dedicated upright facilities

- p-Cure facility at Hadassah Hospital
- operational since 2023
- features synchrotron with fixed beamline
 - synchrotron contains only dipoles ("weak focusing"): compact and simple
- Dedicated in-room upright CT + chair
- Already treated multiple cancer types in patients in clinical pilot studies





Dedicated upright facilities

- Mevion FIT system
 - FDA cleared in September 2025
 - In installation at several sites
- Features very compact cylcotron combined with LEO chair & CT
 - cyclotron is "self-shielded", should fit in an existing photon linac bunker
 - isocentric CT





source: https://mevion.com

Summary upright facilities



- Two very different designs are on the market
- Cyclotron with direct beam + LEO positioner + iso-centric CT
- Synchrotron + robot-mounted chair + in-room CT
- Multiple installations are ongoing
- This is very good for the field!





Some numbers to follow...



Properties of protons and carbon



- protons: an elementary particle, mass $m_0=1{,}673\ 10^{-27}kg$ and charge q = $1.602\ 10^{-19}C$
 - this is ~1800 times heavier than an electron, with same (negative) q
- carbon is an element with 6 neutrons, protons, and electrons
 - for acceleration, has to be ionized by removing up to 6 electrons
 - fully ionized carbon with 6 removed electrons: 12C6+
- a 12C6+ ion has a mass $m_0=1{,}993\ 10^{-26}kg$ and charge ${\bf q}=9.613\ 10^{-19}C$
- These numbers stink so we normalize to unit atomic mass u and unit elementary charge e

Nicer units



- unit elementary charge $e = 1.602 \ 10^{-19} C$
- unit atomic mass $\, {\rm u} = 1.661 \, 10^{-27} kg \,$ conveniently defined as 1/12th of the mass of a carbon atom

- and a unit for energy: electron volts [eV] (remember E=qU)
 - energy of an object with charge ${\bf e}$ at an electric potential of ${\it U}={\bf 1}{\it V}$ ${\rm [eV]}=1.602~10^{-19}{\it J}$
 - using $E=mc^2$, mass can be stated as $m=E/c^2$ in units of [eV]/c²

Properties of protons and carbon



- protons: an elementary particle, mass $m_0=938.3\ MeV/c^2$ and charge ${\rm Q}=1{\rm [e]}$
 - this is ~1800 times heavier than an electron, with same (negative) q
 - the mass of a proton is $m_0=1.007u$ individual protons are heavier due to the mass defect!
- carbon is an element with 6 neutrons, protons, and electrons
 - for acceleration, has to be ionized by removing up to 6 electrons
 - fully ionized carbon with 6 removed electrons: 12C6+
- a 12C6+ ion has a mass $m_0=12u=12~x~931.5~MeV/c^2$ and charge $\rm Q=6~[e]$

Magnetic rigidity of clinical beams



$$B\rho = \frac{m}{q}v = \frac{A[u]}{Q[e]}v$$

Particle	mass at rest m ₀ [MeV/c²]	Mass [u]: Charge [e]	Clinical energies [MeV]	Max Magnetic rigidity [Tm]
Electrons	0.511	0.00055	6 25	0.085
Protons	938.3	1.007	50 230	2.27
Carbon 12C6+	12 x 931.5	2	12 x (80 430)	6.62