

Performance of Prototypes for the PANDA Barrel EMC

Daniel A. Bremer, Tobias Eißner, Peter Drexler, Markus
Moritz, Rainer W. Novotny, C. Rosenbaum for the PANDA
Collaboration

II. Physikalisches Institut, JLU Gießen

FAIR 2014, October 13th – 17th



Bundesministerium
für Bildung
und Forschung

Outline

The PANDA
Detector

PROTO60

Results

Outlook


The PANDA Detector

The Barrel Prototype PROTO60

Test Beam Results

Outlook

Electromagnetic Calorimetry at PANDA





- ▶ Detection of electromagnetic probes with high resolution within a large dynamic range
 - ▶ $10 \text{ MeV} < E < 15 \text{ GeV}$
- ▶ High rate capability
- ▶ Compactness
- ▶ Sufficient radiation hardness
- ▶ Timing information for triggerless DAQ concept

Electromagnetic Calorimetry at PANDA


Daniel
Bremer

The PANDA
Detector
PANDA EMC
PROTO60
Results
Outlook




- ▶ Lead tungstate (PWO-II)
 - ▶ Length 200 mm $\approx 22 X_0$
 - ▶ Operating temperature -25°C


Electromagnetic Calorimetry at PANDA



- ▶ Pb-Sci sampling calorimeter (Shashlyk type)
- ▶ PMT + Sampling-ADC readout





The Barrel EMC Prototype PROTO60



- ▶ Operating temperature -25°C
- ▶ Readout by Sampling ADCs (SIS3302, 50 MHz)
- ▶ Energy and time determined offline with feature extraction algorithms

Setup at CERN-SPS



- ▶ Secondary 15 GeV/c positron beam, 2 positions
- ▶ Calibration via 150 GeV/c muon beam, cosmic MIPs
- ▶ Reference position provided by Si-strip detectors

PROTO60 Scheme

Daniel
Bremer


The PANDA
Detector

PROTO60


Setup

Results

Outlook




Signal Treatment




- ▶ Later on:
Implementation to
FPGAs
- ▶ Sampling rate 50 MHz
⇒ 20 ns Binning

Signal Treatment




- ▶ Later on:
Implementation to
FPGAs
- ▶ Sampling rate 50 MHz
⇒ 20 ns Binning

Raw



MWD




Signal Treatment

The PANDA
Detector

PROTO60
Setup

Results


Outlook




▶ Later on:
Implementation to
FPGAs

▶ Sampling rate 50 MHz
⇒ 20 ns Binning


Raw




MWD +
MA



MWD




MWD +
CFT





SADC Channel number

SADC Channel number



Energy Resolution




- ▶ Overall energy resolution
 - ▶ 15 GeV positron beam
 - ▶ Extracted via fit of Novosibirsk function
- ▶ Position dependence
 - ▶ Restricted beam size $\varnothing 4$ mm



Time Information



- ▶ Time extracted with CFD algorithm
- ▶ Time difference of two neighbouring crystals
- ▶ Energy ratio within 10%
- ▶ Single crystal time resolution:
 $\sigma_t = \frac{\sigma(\Delta t)}{\sqrt{2}} = (441 \pm 15) \text{ ps}$
at 5.5 GeV dep. energy



Position Reconstruction

- ▶ Centre of gravity algorithm


$$x_{cg} = \frac{\sum_i w_i x_i}{\sum_i w_i}$$

- ▶ Logarithmic weighting:

$$w_i = \max \left\{ 0, W_0 + \ln \left(\frac{E_i}{E_T} \right) \right\}$$

$$E_T = \sum_i E_i$$

- ▶ Optimal W_0 requires no further corrections



Position Resolution

The PANDA
Detector

PROTO60

Results


Energy

Time



Position

Summary

Outlook



- ▶ Overall position resolution
 - ▶ $\sigma_x \approx \sigma_y \approx 1.1$ mm
- ▶ Position dependent response
 - ▶ Evaluation in slices of 0.5 mm along x_{S_i}



Summary

- ▶ The PROTO60 was tested at CERN-SPS at the highest design energy of 15 GeV
 - ▶ Obtained Results:
 - ▶ Overall energy Resolution $\frac{\sigma}{E} = 1.7\%$
 - ▶ For restricted beam size ($\varnothing = 4\text{ mm}$) improved resolution up to 1.4%
 - ▶ Time Resolution amounts to $(441 \pm 15)\text{ ps}$ at 5.5 GeV deposited energy
 - ▶ Overall position resolution of $\sigma_x \approx \sigma_y \approx 1.1\text{ mm}$
- ⇒ Set of test beamtimes covering design energy range complete


Outlook - PROTO120

The PANDA
Detector

PROTO60

Results

Outlook
PROTO120



Prototype Comparison

PROTO60

- ▶ 60 crystals, single geometry (type 6)
 - ▶ Single LAAPD, 1 cm² (quadratic)
 - ▶ Low-noise low-power charge preamplifier (LNP)
 - ▶ Trick: different gains in inner and outer ring
- ⇒ Required resolution parameters achieved

PROTO120

- ▶ 120 crystals, three geometries (type 1, 2 and 3)
- ▶ 2 LAAPDs per crystal (rectangular)
- ▶ APFEL 1.4 with final dynamic range
- ▶ close to final mechanics

PROTO120 - First Beam Test at MAMI


Daniel
Bremer

The PANDA
Detector

PROTO60

Results

Outlook
PROTO120



- ▶ Test setup at MAMI with sub-matrix of nine crystals
- ▶ Analysis currently ongoing

Thank you for listening!

