

MR-ToF precision mass measurements of short-lived nuclides

Robert Wolf

– Max-Planck-Institut für Kernphysik –
for the ISOLTRAP Collaboration

NAVI Annual Meeting
17.12.2013, GSI Darmstadt



- *Mass measurements of short-lived nuclides*
- *MR-ToF mass spectrometry*
- *MR-ToF-MS at ISOLTRAP for ...*
- *... isobar separation: ^{82}Zn and Penning-trap stacking*
- *... ion-beam analysis: target/ion source development*
- *... precision mass measurements: n -rich Ca isotopes*
- *Summary*



Mass measurements of short-lived nuclides

$$m(Z, N) = Z \cdot m_p + Z \cdot m_e + N \cdot m_n - BE/c^2$$


→ input for nuclear physics and astrophysics models and calculations

Challenges for experiments on short-lived nuclides:

- half-life: <1s
- production: <100 ions/s
- contaminants: long lived, high yield


Aims and requirements:

- high mass resolving power
- high contamination suppression
- fast measurement cycle
- low measurement uncertainty




Direct mass measurements techniques – the well-known


Penning trap



Storage ring




Time-of-flight- $B\beta$





Multi-Reflection Time-of-Flight Mass Spectrometry

1950s: ToF-MS




$t \sim 0.01 \text{ ms}$, $R \sim 500$

ion flight time:

$$t \propto \sqrt{\frac{m}{q}}$$


mass resolving power:

$$R = \frac{m}{\Delta m} = \frac{t}{2\Delta t}$$






Multi-Reflection Time-of-Flight Mass Spectrometry



$t \sim 0.01 \text{ ms}, R \sim 500$




ion flight time:

$$t \propto \sqrt{\frac{m}{q}}$$


mass resolving power:

$$R = \frac{m}{\Delta m} = \frac{t}{2\Delta t}$$





Multi-Reflection Time-of-Flight Mass Spectrometry




ion flight time:

$$t \propto \sqrt{\frac{m}{q}}$$


mass resolving power:

$$R = \frac{m}{\Delta m} = \frac{t}{2\Delta t}$$





Multi-Reflection Time-of-Flight Mass Spectrometry




ion flight time:

$$t \propto \sqrt{\frac{m}{q}}$$

mass resolving power:

$$R = \frac{m}{\Delta m} = \frac{t}{2\Delta t}$$




Penning traps
decay stations
...



Multi-Reflection Time-of-Flight Mass Spectrometry for short-lived nuclides

Systems in operation




1: Plaß *et al.*, IJMS 349-350, 134 (2013); NIM B 317, 457 (2013); NIM B 266, 4560 (2008); EPJ ST 150, 367 (2007)

2: Schury *et al.*, NIM B 317, 537 (2013); EPJ A 42, 343 (2009); Ito *et al.*, PRC 88, 011306R (2013); Ishida *et al.*, NIM B 241, 983 (2005)

Multi-Reflection Time-of-Flight Mass Spectrometry for short-lived nuclides

Systems in operation and under construction/development




1: Plaß *et al.*, IJMS 349-350, 134 (2013); NIM B 317, 457 (2013); NIM B 266, 4560 (2008); EPJ ST 150, 367 (2007)

2: Schury *et al.*, NIM B 317, 537 (2013); EPJ A 42, 343 (2009); Ito *et al.*, PRC 88, 011306R (2013); Ishida *et al.*, NIM B 241, 983 (2005)

3: Piechaczek *et al.*, NIM B 266, 4510 (2008)

Multi-Reflection Time-of-Flight Mass Spectrometry for short-lived nuclides

Systems in operation and under construction/development



1: Plaß *et al.*, IJMS 349-350, 134 (2013); NIM B 317, 457 (2013); NIM B 266, 4560 (2008); EPJ ST 150, 367 (2007)

2: Schury *et al.*, NIM B 317, 537 (2013); EPJ A 42, 343 (2009); Ito *et al.*, PRC 88, 011306R (2013); Ishida *et al.*, NIM B 241, 983 (2005)


3: Piechaczek *et al.*, NIM B 266, 4510 (2008)



ISOLDE – Isotope Separator On-Line @ CERN

MAX-PLANCK-INSTITUT FÜR KERNPHYSIK



spallation, fission, fragmentation of
UC targets




resonance laser ionization



proton pulses
1.4GeV, 2mA



magnetic
separators
 $\frac{m}{\Delta m} \approx 5000$




Isotope
production:

10^{13}
ions
per
second
 10^{-2}

https://oraweb.cern.ch/pls/isolde/nucl_chart.nuclear_chart




The ISOLTRAP mass spectrometer for short-lived nuclides



The ISOLTRAP mass spectrometer for short-lived nuclides

4 ion traps for ion accumulation, selection, cooling and mass measurement

- 1.) Linear segmented Paul trap (RFQCB)**
ion accumulation and bunching
- 2.) MR-ToF-MS + Bradbury-Nielsen gate**
mass analysis, separation, selection and **mass measurement**
- 3.) Preparation Penning trap**
ion accumulation, cooling and **mass selection**
- 4.) Precision Penning trap**
mass measurement






The MR-ToF-MS at ISOLTRAP

4 ion traps for ion accumulation, selection, cooling and mass measurement


MR-ToF-MS¹:

- mean kinetic energy $E_{\text{kin}} = 2.1 \text{ keV}$
- ToF separation due to different m/q



RFQCB:

- $\Delta t \approx 60 \text{ ns}$
- $\Delta E_{\text{kin}}/E_{\text{kin}} \approx 3\%$



The MR-ToF-MS at ISOLTRAP

MR-ToF-MS

mass resolving power (FWHM)

$m/\Delta m = 100000$ at 12ms

$m/\Delta m = 200000$ at 30ms


transmission

≈50% at 30ms

ion capacity for multiple species

≈1000 per cycle

≈100000 per second





Wolf et al., NIM A 686, 82 (2012); 1: Bradbury & Nielsen, Phys. Rev. 49, 388 (1936); 2: Plass et al., NIM B 266, 4560 (2008)

Bradbury-Nielsen gate (BNG)^{1,2}

contamination suppression

1:10000





Purification of isobaric masses








Isobar separation: n-rich Zn isotopes

MAX-PLANCK-INSTITUT FÜR KERNPHYSIK

MR-ToF spectrum



$N=50$



80Zn:
 $1000/\text{s}, T_{1/2}=0.55\text{s}$
80Rb:
 $10000/\text{s}, T_{1/2}=33.4\text{s}$

Wolf et al., PRL 110, 041101 (2013); NIM A 686, 82 (2012)




Isobar separation: n-rich Zn isotopes



$N=50$

- 0.1s ... 3s
- 3s ... 2m
- 2m ... 1h
- 1d ... 1y
- >1Gy


MR-ToF spectrum



ToF-ICR mass measurement



80Zn:
1000/s, $T_{1/2}=0.55\text{s}$
80Rb:
10000/s, $T_{1/2}=33.4\text{s}$








Isobar separation: n-rich Zn isotopes

$N=50$

0.1s ... 3s
3s ... 2m
2m ... 1h
1d ... 1y
>1Gy



MR-ToF separation in only 2.5ms,
<25ms for complete preparation



82Zn:
200/s, $T_{1/2}=228\text{ms}$
82Rb:
6000/s, $T_{1/2}=1.3\text{min}$

1: Baym, Pethick & Sutherland, ApJ 170, 299 (1971); 2: Pearson *et al.*, PRC 83, 065810 (2011); Wolf *et al.*, PRL 110, 041101 (2013)



$N=50$

0.1s ... 3s
3s ... 2m
2m ... 1h
1d ... 1y
>1Gy

Isobar separation: n-rich Zn isotopes

BPS-model¹ of

neutron-star outer crust:


- masses of n-rich nuclei
- ^{82}Zn predicted by some mass models²
- ➔ first mass measurement:
 ^{82}Zn is not part of the outer crust

ISOLDE techniques:

- resonant laser ionization
- neutron converter
- quartz transfer line

^{82}Zn :
200/s, $T_{1/2}=228\text{ms}$


^{82}Rb :
6000/s, $T_{1/2}=1.3\text{min}$



200ms

15ms

2.5ms




MR-ToF separation in only 2.5ms,
<25ms for complete preparation

1: Baym, Pethick & Sutherland, ApJ 170, 299 (1971); 2: Pearson *et al.*, PRC 83, 065810 (2011); Wolf *et al.*, PRL 110, 041101 (2013)



Isobar separation: preparation Penning trap stacking

- Space-charge effects lead to peak coalescence – limited number of ions
- Multiple MR-ToF-MS separation cycles within one experiment cycle
- Accumulation of isobarically purified bunches in the preparation Penning trap



Isobar separation: preparation Penning trap stacking


- Space-charge effects lead to peak coalescence – limited number of ions
- Multiple MR-ToF-MS separation cycles within one experiment cycle
- Accumulation of isobarically purified bunches in the preparation Penning trap

20Hz repetition rate:

- x40 count rate
- x3.6 meas. time
- **x11 faster meas.**

measurement time

≈1 day → 2 hours!



^{179}Lu :
 $3 \times 10^5/\text{s}$, $T_{1/2}=4.6\text{h}$
 ^{179}DyO :
 $4.5 \times 10^8/\text{s}$, stable




Ion beam analysis



MR-ToF ion-beam analysis

Ion-beam composition analysis



- direct feedback for target/line optimization
- sampling of release curve possible
- single ion sensitivity to detect lowest yields
- no upper limit on half-life as with decay station
- not hindered by decay branching ratio



MR-ToF ion-beam analysis

MR-ToF analyzer to investigate resonant laser ionization of nuclides far from stability

- fast, sensitive tool to improve ionization eff.
- high dynamic range: 1-10e5 counts/s
- counts free from background contamination
- not limited by decay branching ratio
- help to provide isomerically pure beams







Precision mass measurements






MR-ToF mass spectrometer: n-rich Ca isotopes





Wienholtz et al., Nature 498, 346 (2013)



- TITAN/TRIUMF measurements of $^{51,52}\text{Ca}^{(1)}$
- ISOLTRAP measurement agrees with TITAN
- $^{51,52}\text{Ca}$ measured with Penning trap
- $^{52,53,54}\text{Ca}$ measured with MR-ToF-MS
- MR-TOF-MS agrees with PT measurements
- sub-ppm statistical mass uncertainty




1: Gallant *et al.*, PRL 109, 032506 (2012)



Wienholtz *et al.*, Nature 498, 346 (2013)


MR-ToF mass spectrometer: n-rich Ca isotopes

A=53: measurement cycle $\approx 4\text{ms}$



6413 counts/12.6h \rightarrow 9 counts/minute

A=54: measurement cycle $\approx 6\text{ms}$




2314 counts/18.2h \rightarrow 2 counts/minute

statistical uncertainty $\approx 45\text{keV}$ $\rightarrow \delta m/m \approx 9 \times 10^{-7}$

Z=20 40 41 42 43 44 45 46 47 48 49 50 51 52 53 54 55

MR-ToF mass spectrometer: n-rich Ca isotopes


$$S_{2n} = BE(N, Z) - BE(N - 2, Z)$$



Z=20 40 41 42 43 44 45 46 47 48 49 50 51 52 53 54 55

MR-ToF mass spectrometer: n-rich Ca isotopes

$$S_{2n} = BE(N, Z) - BE(N - 2, Z)$$






Z=20 40 41 42 43 44 45 46 47 48 49 50 51 52 53 54 55

MR-ToF mass spectrometer: n-rich Ca isotopes


a

Two-neutron separation energy



c

Two-neutron shell gap




$N=32$ shell closure

excellent agreement with predictions from
microscopic valence-shell calculations
with three-nucleon forces (NN+3N)
based on chiral effective field theory


MR-ToF mass spectrometer at GSI and RIKEN

- first direct mass measurements of $^{211}\text{Rn}^+$, $^{213}\text{Rn}^+$, $^{211}\text{Po}^+$ at GSI
- half-life $^{213}\text{Rn}^+$: 19.5ms
- measurement time 12ms



Plaß *et al.*, NIM B 317, 457 (2013)

- mass measurements of $^8\text{Li}^+$ at RIKEN
- mass resolving power $R \approx 200000$
- measurement time 8ms
- relative mass uncertainty 6.6×10^{-7}




Ito *et al.*, PRC 88, 011306R (2013)

Summary

Precision and fast measurement cycle makes the MR-TOF-MS a promising approach for MS on short-lived isotopes with low production

MR-ToF mass purification for Penning trap mass spectrometry

- $R=200000$ after 30ms
- contam. suppression 10000:1
- first mass measurement of ^{82}Zn



Target/ion source development and optimization facilitated

MR-ToF mass measurements on very short-lived isotopes with very small production rates

- successful mass measurements of $^{52-54}\text{Ca}$, sub-ppm uncertainty
- in addition: $^{52-53}\text{K}$, half-life 30ms, publication in preparation

Other groups and activities (for short-lived nuclei):

FRS ioncatcher (GSI, Gießen), SlowRI (RIKEN), ORNL (Oak Ridge), CARIBU (Argonne), IGISOL (Jyväskylä), GANIL (Caen), TITAN/TRIUMF (Vancouver), RISP (S. Korea)

SPONSORED BY THE

Thanks to: the ISOLTRAP collaboration

P. Ascher, D. Atanasov, D. Beck, K. Blaum, Ch. Böhm, Ch. Borgmann, M. Breitenfeldt, R. B. Cakirli, T. E. Cocolios, S. Eliseev, T. Eronen, D. Fink, S. George, F. Herfurth, A. Herlert, S. Kreim, M. Kowalska, Yu. Litvinov, D. Lunney, V. Manea, E. Minaya-Ramirez, S. Naimi, D. Neidherr, M. Rosenbusch, L. Schweikhard, J. Stanja, F. Wienholtz, R.W., K. Zuber

the ISOLDE collaboration, CERN



<http://isoltrap.web.cern.ch>



Federal Ministry
of Education
and Research
Grants No.:
05P12HGCI1
05P12HGFNE

