

I. Physikalisches <u>JUSTUS-I</u> Institut GIE





MR-TOF-MS for Nuclear Physics Facilitate a New Class of Experiments in Analytical Mass Spectrometry

T. Dickel^{1,2}, J. Lang¹, W. Lippert¹, W.R. Plaß^{1,2}, S. Ayet San Andrés^{1,2}, J. Bergmann¹, H. Geissel^{1,2}, A. Pikhtelev³, C. Scheidenberger^{1,2}, M.I. Yavor⁴

¹ II. Physikalisches Institut, Justus-Liebig-Universität Gießen, Germany
² GSI Helmholtzzentrum für Schwerionenforschung, Darmstadt, Germany
³ Institute for Energy Problems of Chemical Physics, RAS Chernogolovka, Russia
⁴ Institute for Analytical Instrumentation, RAS, St. Petersburg, Russia

Overview

- Motivation
- MR-TOF-MS for Analytical Mass Spectrometry
- Highly Resolving Tandem MS
- Applications
- Summary

Mobile High-Performance Mass Spectrometry

High-Performance Mass Spectrometers

Very high resolving power (m/ Δ m > 10⁵)

- Resolution of isobars at low/medium mass
- Isotopic resolution at high mass
- Very high mass accuracy ($\delta m/m < 10^{-6}$)
- Determination of composition and structure

Mobile Mass Spectrometers

- In-situ measurements
- Prompt results
- Investigation of time-dependent processes
- No sample storage required



A. Marshall et al., 9.4 T FT-ICR-MS



R.G. Cooks et al., Mini 11 RIT-MS

There is no mobile high performance MS yet Goal: develop mobile high performance MS

MR-TOF-MS

Different fields of application

Synergy

Fast, sensitive, broadband, non-scanning, high resolving power

MR-TOF-MS for nuclear physics experiments



W.R. Plaß et al., NIM B 266 (2008) 4560 W.R. Plaß et al., Int. J. Mass Spectrom. 394 (2013) 134 T. Dickel et al., NIM A 777 (2015) 172

> J. Ebert, Mo. 12:40 C. Hornung, Tu. 18:30 S. Ayet, Tu 18:30

No magnet, moderate vacuum requirements, size reduction is straight forward

MR-TOF-MS for in-situ analytical mass spectrometry





Water monitoring at 'hot spots'

Mobile MR-TOF-MS for Analytical Mass Spectrometry



Electrospray Ionisation

Soft ionisation technique for large molecules, e.g. proteins

Multistep process:

- 1. Analyte in solvent
- 2. Spray solvent towards the inlet of the mass spectrometer
- 3. Charged solvent droplets evaporate and intact charged macromolecule remain



Desorption Electrospray Ionisation

Z. Takats / R.G. Cooks, 2004



Mobile MR-TOF-MS for Analytical Mass Spectrometry

Design goals

- Mobile
- $m/\Delta m > 100,000$
- δm/m < 1 ppm
- Repetition rate ~ kHz
- Compatible with atmospheric ion sources

Rugged Setup

- 3 stacked DN150CF vacuum crosses as recipient
- 30 m³/h rough pump integrated
- Volume ~ 0.8 m³
- Weight ~200kg
- No further infrastructural needs except one power line



T. Dickel et al, *Nim B*, **317**, 779-784 (2013) J. Lang, PhD thesis in preparation, JLU Gießen

Time-of-Flight Analyzer



¹M.I. Yavor, et al., Int. J. Mass Spectrom. (2015) in press, doi:10.1016/j.ijms.2015.01.002

Tandem Mass Spectrometry

Problem:

- Unambiguous identification of bio-molecules despite high mass multiplicity
- Structural information

Solution: Tandem-MS - two or more stages of mass analysis in space or time



Advantages of MS^N with high resolving power in all steps:

- \rightarrow no preseparation (GC, LC, IMS) neccessary
- \rightarrow analysis of more complex mixtures
- \rightarrow higher confidence in identification

Ideal for In-situ analysis

Mass Selective Retrapping

1. Mass selective ion retrapping



- 1. Separate ions in time-of-flight
- 2. Inject ions back in RF-Trap
- 3. Close RF-Trap
 - \rightarrow only ions of interest are recaptured



Collissional-Induced Dissociation

2. Collisional-induced dissociation in the ion trap

What is collisional-induced dissociation (CID)?

An ion/neutral species interaction wherein the projectile ion is dissociated as a result of interaction with a target neutral species.



"Low-energy" CID

Fime scale	
Collision number:	
Collision energy:	

~ ms 10-100 1-200eV

CID by resonant excitation in the ion trap

High Resolving Tandem Mass Spectrometry

1-3. Ion re-trapping \rightarrow CID \rightarrow mass analysis



Application I: Water-Monitoring at "Hot Spots"

Check and monitor wastewater for pharmaceuticals (eg. Carbamazepine) and antibacterial agents (eg. Triclosan)

Two-step process:

1. Enrichment of analyte by thin-film microextraction

2. Sample analysis by high resolution DESI-MS



Application II: Real-time Tissue Recognition

Problem:

• Histological examination of tissue is slow (hours)

Solution:

Rapid Evaporative Ionization Mass Spectrometry (REIMS) Identification of healthy and cancerous tissue in less than a second by the mass spectra





Vapor generated by a electrosurgery is pumped to API of mass spectrometer



Z. Takats, et al., Anal. Chem. 2010, 82, 7343-7350

Summary

- Motivation: High-resolution in-situ mass spectrometry Method: Mobile MR-TOF-MS
- Device has been developed and commissioned
 - $m/\Delta m = 4*10^{5}$
 - Atmospheric Pressure Inlet with various ionization techniques
- Highly resolving Tandem MS Mass selective Retrapping + CID → MS^N with high resolution >10⁴ in all steps
- in-situ applications planned eg. Wastewater monitoring
- Next generation device: portable
 - Replace standard components by customized components
 - \rightarrow further saving in size, weight, power consumption, ...
- Future applications:
 - Space missions





Acknowledgments

IONAS Group at JLU Gießen

S. Ayet, J. Bergmann, U. Czok, T. Dickel, M. Diwisch, J. Ebert, H. Geissel, F. Greiner, E. Haettner, S. Heinz, C. Hornung, C. Jesch, R. Knöbel, N. Kuzminchuk, J. Lang, M. Petrick, A. Pikhtelev, W.R. Plaß, M.P. Reiter, A.-K. Rink, C. Scheidenberger, B. Sun, W. Lippert, M. Werner, M. Yavor





Financial support was provided by the BMBF under contracts No. 05P12RGFN8, by the state of Hesse (LOEWE-Schwerpunkt AmbiProbe and LOEWE Center HICforFAIR).

JUSTUS-LIEBIG-UNIVERSITÄT GIESSEN