

Development of a rapid solvent extraction apparatus for aqueous chemistry of the heaviest elements

Yukiko Komori for a RIKEN – Niigata Univ. – JAEA – Univ. Tsukuba – Tohoku Univ. – Univ. Oslo collaboration (The SHE aqueous chemistry collaboration at GARIS)

Introduction: Aqueous chemistry of SHEs

Chemistry of SHEs

104	105	106	107	108	112	113	114
Rf	Db	Sg	Bh	Hs	Cn		Fl

Gas: $Z = 104-108, 112-114$

Aqueous: $Z = 104-106$

Nuclide	Half-life	Production rate*
$^{261}\text{Rf}^a$	68 s	420 atoms/h
^{262}Db	34 s	70 atoms/h
$^{265}\text{Sg}^{a,b}$	8.5 s/14.4 s	12 atoms/h
^{266}Bh	10.7 s	1.7 atoms/h

* ^{248}Cm target thickness: 300 $\mu\text{g}/\text{cm}^2$; Beam intensity: 2 μA

Pioneering cation-exchange studies of Sg in HNO_3/HF and HNO_3

Schädel et al., Radiochim. Acta **77**, 149 (1997).; Radiochim. Acta **83**, 163 (1998).

Conventional aqueous chemistry apparatus used for Rf, Db, and Sg

ARCA and AIDA: batch-wise column chromatography apparatuses with Si detectors for α/SF spectrometry


- Decay loss during aerosol collection (~ 30 s)
- Decay loss during α -source preparation (~ 30 s)
- Low detection efficiency: $\text{eff.}(\alpha) = \sim 30\%$
→ $\text{eff.}(\alpha-\alpha) = \sim 9\%$; $\text{eff.}(\alpha-\alpha-\alpha) = \sim 3\%$
- A huge amount of background radioactivities of by-products

RIKEN GARIS gas-jet system

Requirements for aqueous chemistry studies of Sg and the heavier SHEs:

- **Continuous and rapid chemical separation**
- **Rapid and efficient α /SF detection under low-background condition**


GARIS gas-jet system is ready for SHE chemistry at RIKEN:



- By-products can be removed almost completely.
- Liquid scintillation (LS) detectors with a high detection efficiency ($\sim 100\%$) will become available for aqueous chemistry of SHEs.

Purpose of this study

Development of a continuous and rapid solvent extraction apparatus coupled to the GARIS gas-jet system for aqueous chemistry of the heaviest SHEs



Continuous dissolution (MDG), solvent extraction (FSE), and radiation detection with a flow LS detector

- **Rapid chemical separation and α -source preparation**
→ **Minimum decay loss**
- **High-detection efficiency (~100%) for α - α and α -SF correlations**

Feasibility of aqueous chemistry of Sg and Bh

Production and decay studies of $^{265}\text{Sg}^{a,b}$ ($T_{1/2} = 8.5 \text{ s}, 14.4 \text{ s}$) and ^{266}Bh (10.7 s):

$^{248}\text{Cm}(^{22}\text{Ne}, 5n)^{265}\text{Sg}^{a,b}$ Haba et al., Phys. Rev. C **85**, 024611 (2012).

$^{248}\text{Cm}(^{23}\text{Na}, 5n)^{266}\text{Bh}$ Haba et al., TASCA15 contribution.

Continuous solvent extraction and LS detection (Present apparatus)

Nuclide	$T_{1/2}$ [s]	σ [pb]	Target [$\mu\text{g}/\text{cm}^2$]	Beam [p μA]	Cool. T. [s]	Chem. Y. [%]	Detec. eff.* [%]	Event rate [/d]
$^{265}\text{Sg}^a$	8.5	180	300	4	10	50	100	3.2
$^{265}\text{Sg}^b$	14.4	200	300	4	10	50	100	5.3
^{266}Bh	10.7	55	300	4	10	50	100	1.2

Batch-wise chemical separation (e.g. ARCA and AIDA)

Nuclide	$T_{1/2}$ [s]	σ [pb]	Target [$\mu\text{g}/\text{cm}^2$]	Beam [p μA]	Coll. T. [s]	Cool. T. [s]	Chem. Y. [%]	Detec.* eff. [%]	Event rate [/d]
$^{265}\text{Sg}^a$	8.5	180	300	4	30	30	50	9	0.02
$^{265}\text{Sg}^b$	14.4	200	300	4	30	30	50	9	0.1
^{266}Bh	10.7	55	300	4	30	30	50	9	0.01

* Efficiencies for α - α correlations.



This study

- **Development of Membrane DeGasser (MDG) and Flow Solvent Extractor (FSE)**
- **Performance evaluation of MDG and FSE**
- **Online solvent extraction of Tc and Re with MDG-FSE**

Development (1): MDG

Univ. Oslo/JAEA Membrane DeGasser (MDG)

Ooe et al., J. Radioanal. Nucl. Chem. **303**, 1317 (2015).



Dissolution efficiency of $^{91\text{m}}\text{Mo}$ ($T_{1/2} = 65$ s):

- > 80% at high flow rates of 6–24 mL/min
- decreases with a decrease of the aq. flow rate.
50–60% at a lower flow rate of 1 mL/min



Development (1): MDG

RIKEN Membrane DeGasser (RIKEN-MDG)

A new MDG was fabricated by modifying Univ. Oslo/JAEA-MDG to dissolve shorter-lived nuclides with high efficiencies at a low flow rate of ~ 1 mL/min.

Major modifications:

- Dead volume: $\sim 90 \mu\text{L} \rightarrow \sim 23 \mu\text{L}$
- Static mixer \rightarrow Simple T-connector







Development (2): FSE

Flow Solvent Extractor (FSE)




Karlberg and Thelander, Anal. Chim. Acta **98**, 1 (1978).

Bergamin et al., Anal. Chim. Acta **101**, 9 (1978).


Motomizu and Oshima, Analyst **112**, 295 (1987).



Experimental (1): Performance evaluation of MDG





Nuclear reactions:



Experimental (2): Performance evaluation of FSE

- Production of long-lived and no-carrier-added radiotracers at RIKEN AVF: $^{nat}\text{Mo}(d,xn)^{95m}\text{Tc}$ ($T_{1/2} = 61$ d) and $^{nat}\text{W}(d,xn)^{183}\text{Re}$ ($T_{1/2} = 70$ d)
- Extraction with FSE: HNO_3 -Tri-*n*-octylamine (TOA) / toluene




→ Determination of distribution ratio, $D = [A]_{\text{org.}}/[A]_{\text{aq.}}$; A: radioactivities


	<i>D</i> vs. Capillary length	<i>D</i> vs. [TOA]
Aq. phase	0.1, 1 M HNO_3 + ^{95m}Tc , ^{183}Re	1 M HNO_3 + ^{95m}Tc , ^{183}Re
Org. phase	0.01 M TOA / toluene	0.01, 0.05, 0.1 M TOA / toluene
Capillary length	5, 10, 20, 30, 40, 50, (60), 100 cm	100 cm

→ Comparison with D in equilibrium in the batch extraction (30-min shaking)

Experimental (3): Online solvent extraction of Tc and Re with MDG-FSE




Nuclear reactions:



- **FSE ext. (1): D vs. Capillary length, L**
 $L = 5, 10, 20, 30, 40, 50, 70,$ and $100\ \text{cm}$
- **FSE ext. (2): D vs. [TOA]**
 $[\text{TOA}] = 0.005, 0.01, 0.05,$ and $0.1\ \text{M}$
- **Batch ext. (3-min shaking)**

Results and discussion (1): Performance of MDG

γ -ray spectrum of an aqueous solution dissolved with MDG



Dissolution efficiency with MDG

Nuclide	$T_{1/2}$	Dissolution eff.*
^{90m}Nb	18.8 s	$56 \pm 2\%$
^{90g}Nb	14.6 h	$88 \pm 6\%$
^{178a}Ta	2.36 h	$82 \pm 7\%$


* He gas: 1.5 L/min; 1 M HF: 1 mL/min

- The dissolution efficiency of $\sim 60\%$ was obtained with RIKEN-MDG for the short-lived ^{90m}Nb even at a low aq. flow rate of 1 mL/min.
- Reduction of chemicals and radioactive wastes
- Reduction of quenching effects and increase of energy resolution in α/SF -spectrometry with a LS detector.

Results and discussion (2): Performance of FSE


D vs. Capillary length

[TOA] = 0.01 M



D vs. TOA concentration

Capillary length = 100 cm




- **Extraction equilibrium is attained with the 40-cm capillary.**
Time required for solutions to pass through the 40-cm capillary: ~2.4 s
- **D values with FSE consistent with those by the batch method.**
- **FSE is applicable to determine D values in the wide D range:**
D = ~0.1 – ~20.

Results and discussion (3): Online solvent extraction of Tc and Re with MDG-FSE


D vs. Capillary length

[TOA] = 0.01 M



D vs. TOA concentration

Capillary length = 100 cm



- Discrepancies in D values between FSE and the batch extractions were found for $^{92,94}\text{Tc}$ at $[\text{TOA}] > 0.05 \text{ M}$.

- Online solvent extraction of Tc and Re was successfully performed with stable and high chemical yields: $92 \pm 3\%$ (^{181}Re) during the 6-h beam time

Summary

- We have developed a new rapid chemistry apparatus which consists of **MDG** and **FSE** for the aqueous chemistry studies of Sg and Bh at GARIS.
- Online solvent extraction of Tc and Re was successfully performed with **MDG-FSE** in HNO₃-TOA/toluene.
 - Rapid extraction equilibrium: ~2.4 s (40-cm capillary)
 - Wide applicable D range: $D = \sim 0.1 - \sim 20$
 - High chemical yield: > 90% (¹⁸¹Re)
 - Stable running: > 6 h
 - Low flow rate: 1 mL/min
- **A flow liquid scintillation detector** will be developed by referring to the knowhow from SISAK.
- Interesting chemistry systems for Sg and Bh are under study using radiotracers of their homologues.

Collaborators for the aqueous chemistry at GARIS

Nishina Center for Accelerator-Based Science, RIKEN

H. Haba, S. Yanou, and K. Watanabe

Niigata Univ.

K. Ooe, M. Murakami, D., Sato, and R. Motoyama

Advanced Science Research Center, JAEA

A. Toyoshima and A. Mitsukai

ELPH, Tohoku Univ.

H. Kikunaga

Univ. Tsukuba

A. Sakaguchi and J. Inagaki

Univ. Oslo

J. P. Omtvedt