Second Workshop on Nuclear Spectroscopy Instrumentation Network of ENSAR2

NUSPIN 2017

GSI Darmstadt 28.06.2017

New expertise in γ - γ fast timing using arrays of LaBr₃(Ce) detectors

SPONSORED BY THE



Jean-Marc Régis

Institute of Nuclear Physics University of Cologne, Germany





Content:

The Generalized Centroid Difference Method

to analyze **γ–γ time-difference spectra** from a large **fast-timing array**

The EXILL&FATIMA campaign 2013

at the Institut Laue-Langevin



Results of **Germanium-gated** γ**–**γ **fast timing**

of excited states in **fission fragments**


How to correct for background contributions to the time spectra?


Basis: The picosecond sensitive Centroid Shift Method

The prompt response function (PRF) obtained for τ < 1 ps.

The centroid or center of gravity is the first moment of a time distribution D(t):

$$C^{D} = \langle t \rangle = \frac{\int t D(t) dt}{\int D(t) dt}$$






Asuuming no background contributions, the **centroid** of the "delayed" time spectrum D(t) is **shifted** from the centroid of the convoluted PRF P(t) by the mean lifetime τ of the excited state:

$$\tau = C^{D}(E_{\text{start}}, E_{\text{stop}}) - C^{P}(E_{\text{start}}, E_{\text{stop}})$$

"the time-walk characteristics"

In general: $C^P(E_{\text{start}}, E_{\text{stop}}) = T_0 + T(E_{\text{start}}) + T(E_{\text{stop}})$, where T(E) is the detector time response.

The γ - γ fast-timing technique and the Generalized Centroid Difference method




The **Prompt Response Difference** (PRD) describes the linearly combined γ – γ time walk of the fast-timing array.

J.-M. Régis et al., NIM A **726** (2013)


The standard PRD calibration procedure using the ¹⁵²Eu γ-ray source E [keV] level τ [ps] - nearly no background (no correction needed) 1643.4 - precisely known picosecond lifetimes 3_{1}^{+} - several γ rays connecting one state for different states 1434.0 778.9 367.8 a LaBr₃ gate on 344 keV 779 1123.2 8000 PRD(779,344) = 51(5) psp/b=42(3)FWHM = 405(3) ps930.5 10.5(9)6000 counts/10 ps counts/keV anti-delayed delayed (344,779)(779,344)10.5(6) 755.4 4000-10 p/b=17(2)2000 1089 1299 344.3 46.7(2.5) III II 800 900 1000 1100 1200 1300 0.5 1.5 t [ns] $E_{\gamma}[\text{keV}]$ b LaBr₃ gate on 444 keV $\tau_{Lit} = 1.3(1) \text{ ps}$ $\Delta C = -88(4) \text{ ps}$ ^{152}Gd PRD(444,964) = -91(4) ps100 FWHM = 358(4) psp/b=30(3)E [keV] counts/keV counts/10 ps delayed anti-delayed 1085 (444,964)(964,444) 1579.4 1529.8 ~40 fs

a)



Picosecond sensitive lifetime determination using the <u>mirror-symmetric</u> GCD method.

The PRD calibration procedure using ¹⁵²Eu:



LaBr₃-LaBr₃ coincidences with virtually no background contributions.


no correction and therefore, no systematic error is introduced.

Advantage of the GCD method:

PRD uncertainty within 3**σ**: typically 5-10 ps

The Generalized Centroid Difference (GCD) Method for y-y fast-timing arrays

An analog γ – γ fast-timing array circuitry for trigger-less data acquisition:




$$\{1,2,3,4\} \times \{5,6,7,8\} \in TAC A$$

$$\{1,2,5,6\} \times \{3,4,7,8\} \in TAC B$$

$$\{1,3,5,7\}$$
 x $\{2,4,6,8\}$ \in TAC C

The TAC_{ii} projections:



Off-line alignment of the TAC_{ij} projections using shift_{ij} constants and superposition of the data to perform the GCD method.

Only combinations with i < j are accepted (simplified sorting algorithm). Invalid or multiple combinations are excluded off-line. => 3 TAC and 2 FAN modules for 8 detectors (28 combinations).

Main advantages: no dead-time effects for detector rates of up to 30 kHz and almost no degradation of time resolution (<10 ps).

The EXILL&FATIMA mixed array for prompt γ-ray spectroscopy of cold-neutron induced capture/fission experiments at ILL 2013

Ring of 8 BGO shielded EXOGAM clovers used to provide one or two selective γ-triggers.


16 almost equal LaBr₃(Ce) detectors for γ – γ lifetime measurements.

Collimated cold-neutron beam

Ø1.2 cm

Trigger-less digital data acquisition of 71 digitiser channels

Detector rate: up to 25 kHz Data rate: up to 6.5 MB/s Acquired data: ~ 40 TB




Target position

FATIMA: Inter-detector Compton scattering (cross talk) detected in neighbouring LaBr₃(Ce)detectors.


Off-line elimination (loss of statistics: about 15%). Small Compton suppression of about 10% achieved.

Energy performance of the EXILL & FATIMA spectrometer @ ILL 2013



Timing performance of FATIMA @ ILL 2013

Consited of: 8 cylindrical Ø1.5"x1.5" and 8 cylindrical Ø1.5"x 2" LaBr3(Ce) scintillators




For τ < FWHM:

The statistical centroid uncertainty: $\delta C \sim \frac{FWHM}{2.355\sqrt{r}}$

e.g. FWHM=500 ps and n=1000 counts: $\delta C \sim 7$ ps

Uncertainty of the PRD determination corresponding to a 3σ standard deviation: $\delta PRD = 10$ ps





EXILL&FATIMA 2013:

Identification of fission fragments & investigation on background and contamination

10 days cold-neutron-induced fission of ²³⁵U


fisson yield of 98 Sr: 0.81%





The LaBr3 peak-to-background ratio in this case is smaller than 0.5.

About 65% of the events are related to the background underneath the two FEPs of the γ - γ cascade.




The first 4+ state in ⁹⁸Sr: τ_{lit} = 115(9) ps (β–γ timing, H. Mach et al., Phys. Lett. B 230, 1989)

Assuming only one background contribution, this correction is exact: $t_{cor.} = \frac{\Delta C_{exp.} - \Delta C_{BC}}{P/B}$ Using γ - γ timing: The time responses ΔC_{BC} of the two background components are needed.

As a good approximation:

$$\Delta C_{\text{FEP}} = \Delta C_{\text{exp.}} + [t_{\text{cor.}}(E_{\text{feeder}}) + t_{\text{cor.}}(E_{\text{decay}})]/2$$




 $\Rightarrow \Delta C_{FEP} = 202(12) \text{ ps}, PRD = -41(7) \text{ ps} \Rightarrow \tau = (\Delta C_{FEP} - PRD)/2 = 121(11) \text{ ps}$

The first 2+ state in 94 Sr: $\tau_{lit} = 10(4)$ ps $(\beta-\gamma \text{ timing, H. Mach et al., NPA 523, 1991})$

Experimental observations:

Generally, the shifts of the two background components relative to the total time spectrum are different, sometimes with opposite sign.



counts per 10 ps

10

 $\Delta C_{\text{exp.}} = -31(8) \text{ ps}$

837-1309

anti delaved

1309-837

delayed

Considering the Slope Method for long lifetimes $\tau > 1$ ns

The first 2⁺ state in 98 Sr: $\tau_{lit} = 4.01(12)$ ns $(\beta - \gamma \text{ timing, H. Mach et al., Phys. Lett. B 230, 1989})$

In case of a two component time distribution, the most precise result for the longer decay


220

240

140

E_γ [keV]

120



Conclusions

The **EXILL&FATIMA** campaign was a ground-breaking facility for future γ-γ timing experiments using **large fast-timing arrays** such as: **DESPEC for FAIR@GSI**.

So far: More than 30 lifetimes in more than 12 nuclei have been measured in the range of 10 ps to 4 ns. 8 scientific articles on nuclear structure physics have been published.

Precise and consistent values could be determined for cases with peak-to-background ratio of 0.3 to 2.5.

New GCD method to generate and analyze γ - γ time distributions. Good background-time correction by precise determination of peak-to-background ratios and centroids of background-time distributions.

Experimental observations on the timing behavior of the background are confirmed by recent GEANT4 Monte Carlo simulations.


Thank you very much for your attention

Monte Carlo Simulation of a γ - γ fast-timing experiment using GEANT4


- 2 LaBr3(Ce) detectors face-to-face with distance of 10 cm
- 3 different γ-ray interactions according to their cross sections
- location of each interaction and the corresponding energy release recorded
- intrinsic zero-time response related to full-energy peak (FEP) events
- **TEST**: non-linear time difference between FEP and Compton events (60 Co);


simulation vs. experiment: exact agreement!

Simulation:





M. Dannhoff, IKP Köln, PhD Thesis in work





Monte Carlo Simulation of a γ - γ fast-timing experiment using GEANT4





Component	centroid	relative intensity
Total (Exp.)	402	100%
FEP – FEP	408	47%
FEP – BG BG – FEP BG – BG	392 412 393	30% 12% 11%

First observation using different energies, lifetimes and P/B ratios:

No significant systematical error could be measured using the proposed time correction given on slide 12.

to be continued ...

Generated using a constant $\boldsymbol{E}_{\text{feeder}}$ - $\boldsymbol{E}_{\text{decay}}$ condition.

M. Dannhoff, IKP Köln, PhD Thesis in work